

Safe Job Analysis (SJA)

SJA title: High pressure fluid measurements		Project no: Generic		Project title: Generic	
Project manager: Kjetil Folgerø		Responsible for SJA: Anna Mathea Skar		Participants in the SJA: Anna Mathea Skar, Kjetil Folgerø	
Short work description Use of pressure vessel to test permittivity of non-flammable fluids				Premises: NORCE Bergen	
No.	Task	Possible danger	Possible consequence	Recommended action	
1	Hydrostatic pressure test of measurement cell	<ul style="list-style-type: none">• Deformation of equipment that cannot withstand the pressure• Water leakage	<ul style="list-style-type: none">• Impact damage to personnel or equipment• Damage to electrical equipment and personnel	<ul style="list-style-type: none">• Follow the order of procedure (attached) and checklist (below)• Use mechanical separation between operator and measuring cell• Use safety glasses• Carry out experiments during working hours and under the supervision of a laboratory supervisor who has sufficient training for the equipment• Ensure communication with the project team, laboratory users and those affected nearby before start (and when finished)• Check that all connections to the measuring cell are securely fastened• Consider whether the measuring cell and cables should be permanently installed• Interrupt measurement in the event of a major water leak.• Have fire safety equipment available• Follow a gradual pressure build-up and continuous check of the pressure gauge• Release the pressure before removing the mechanical separator• Do not put your head over the lid when opening the cylinder• Choose measures to avoid exposure to hazards, i.e. ensure that the working environment is suitable, e.g. work room, ventilation, first aid equipment	

Safe Job Analysis (SJA)

2	Initial installation	<ul style="list-style-type: none"> Gas leakage may occur High pressure produces large forces that can destroy equipment 	<ul style="list-style-type: none"> Health damage to personnel by inhaling gas Impact damage to personnel or equipment 	<ul style="list-style-type: none"> Follow procedure and checklist Ensure that the operator has sufficient training/knowledge Follow the emergency shutdown procedure (below) in case of leakage or other unwanted incidents Consider the use of mechanical separation between operator and measuring cell Use safety glasses Ensure communication with the project team, laboratory users and any others affected before start (and when finished) Check that all connections to the measuring cell are securely fastened Make sure that the gas outlet is discharged into a fume hood Consider whether the measuring cell should be placed under extraction Consider whether the measuring cell and cables should be fixed
3	Gas testing	<ul style="list-style-type: none"> Gas leak Deformation of equipment 	<ul style="list-style-type: none"> Health damage to personnel by inhaling gas Impact damage to personnel or equipment 	<ul style="list-style-type: none"> Follow procedure and checklist Ensure that the operator has sufficient training/knowledge Follow the emergency shutdown procedure (below) in case of leakage or other unwanted incidents Use safety glasses Clear sign "gas under pressure" Consider how high a pressure is necessary for the gas being measured. A goal should be to use as low a pressure as possible If the gas leak is out of control, leave the laboratory and notify the room manager. In case of serious personal injury, call 113 and the NORCE emergency number (21 08 01 86)
4	Disconnection	<ul style="list-style-type: none"> Gas leak Deformation of equipment 	<ul style="list-style-type: none"> Health damage to personnel by inhaling gas Impact damage to personnel or equipment 	<ul style="list-style-type: none"> Follow procedure and checklist Use safety glasses Ensure communication with the project team, laboratory users and any others affected when finished Release the pressure and ensure that the pressure has dropped to atmospheric pressure before removing the mechanical separation between operator and measuring cell

Check list:

	Yes	No	Comments
Have you informed the room manager and the person in charge of emergency response for the laboratory?			Inform room manager well in advance of the start of the experiment. Also inform when the experiment is finished
Will the experiment be done under supervision?			Consider supervision at high pressures and if there is a lack of sufficient training and experience with the equipment used
Are the fluid(s) explosive?			Follow a different/create new SJA and always conduct experiments with supervision
Are the fluid(s) health hazardous?			Be careful using gas exhaust hoods during gas emissions. Use a laboratory room with adequate ventilation
Should the fluid(s) be pressurized?			Follow the procedure and be careful with the order of procedure
Is the gas from the measuring cell released into the exhaust before the chamber is opened or the separator is removed?			Check that the pressure sensor has returned to atmospheric pressure

Emergency shutdown procedure:

1. If the gas leak is out of control, leave the laboratory and notify the room manager.
2. In case of serious personal injury, call 113 and the NORCE emergency number: 21 08 01 86